

Data Sheet (DS) 7TM0318N

Issuing Date: 2020-05-18

Description	KOP (non-phospho), κ-Opioid Receptor Antibody	
Format	Purified, Liquid	
Product Type	Rabbit Polyclonal Antibody	
Isotype	Polyclonal IgG	
Quantity	100 μΙ	

Product Details

Applications	This product has been reported to work in the following applications:		
		Dilution	
	Western Blot	1:1000	
	This information is derived from testing within our laboratories and peer-reviewed publications. Please refer to references indicated for further information. For general protocol recommendations, please visit https://Ttmantibodies.com/7tm-antibodies-support/7tm-protocols/ Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.		
Target Species	Human, Mouse, Rat		
Product Form	Purified IgG, liquid		
Antiserum Preparation	Antiserum to κ-Opioid Receptor was raised by repeated immunization of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.		
Immunogen	A synthetic peptide presents carboxyl-terminal tail of human KOP.		
Storage Buffer	Dulbecco's PBS, pH 7.4, with 150 mM NaCl, 0.02% sodium azide		
Specificity	The non-phospho-kappa-opioid receptor antibody is directed against the distal end of the carboxyl-terminal tail of mouse, rat and human KOP. In can be used to detect total MOP receptors in Western blots independent of phosphorylation.		
Guarantee	12 months from date of dispatch		
Storage	Store at -20°C. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody.		
Regulatory	For research purposes only		
Health and Safety Information	Material Safety Data Sheet documentation is available at		

Safety Data Sheet EU Safety Data Sheet US

Related Products

pS356/pT357-KOP (phospho-к-Opioid Receptor Antibody)

pT363-KOP (phospho-к-Opioid Receptor Antibody)

Details of the Supplier of the Data Sheet

Supplier

7TM Antibodies GmbH

Hans-Knöll-Str. 6

07745 Jena – Germany

TEL: ++49 151 20130575

FAX: ++49 3641 241 49 58

Email: info@7tmantibodies.com

Website: 7tmantibodies.com